METTLER TOLEDO

**InPro® 5500 i
CO₂ Sensor**

**Instruction manual
Bedienungsanleitung
Instructions d'utilisation**

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Contents

Inhalt

1	Introduction	5
2	Important notes	6
2.1	Notes on operating instructions	6
2.2	Intended use	6
2.3	Safety instructions	7
2.4	Examples of some typical applications	8
3	Product description	9
3.1	General information	9
3.2	Principle	9
3.2.1	Thermal conductivity measurement of CO ₂ in liquids	9
3.2.2	Determination of CO ₂ concentration	10
3.2.3	Measurement cycles with the InPro 5500 i	10
3.3	Scope of delivery	11
3.4	Equipment features	12
4	Installation	13
4.1	Mounting the sensor	13
4.1.1	Mounting and positioning the sensor	13
4.1.2	Removing the sensor and storage	13
4.2	Purge air supply	14
4.2.1	Requirements of the purge gas supply:	14
4.2.2	Connecting the purge air to the sensor	14
4.3	Connections	15
4.2.1	Connecting the sensor to an RS-485 cable	15
4.2.2	Connecting the cable to the transmitter	15
5	Operation	16
5.1	Start-up	16
5.2	Configuration	16
5.2.1	Sensor detection	16
5.2.2	CO ₂ solubility and temp. factor	16
5.3	Calibration	17
5.3.1	Purpose of calibration	17
5.3.3	Single point calibration in-line, off-line	17
5.3.4	In-line process single-point calibration	18
6	Maintenance	19
6.1	Inspection of the sensor	19
6.1.1	Visual inspection	19
6.1.3	Testing the sensor with a transmitter	19
6.1.4	ISM diagnostics	20
6.1.5	ISM design	21
6.2	Changing the MembraCap	21
6.3	Error handling	25
7	Storage	26
8	Product specification	27
8.1	Certificates	27
8.2	Specifications	28
9	Ordering information	29
9.1	Sensors	29
9.2	Accessories	29
9.3	Spare parts	29
9.4	Recommended transmitters	29

1 Introduction

Thank you for buying the **InPro 5500 i sensor from METTLER TOLEDO.**

The construction of the InPro series employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment.

Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.

The instruction manual must always be stored close at hand, in a place accessible to all people working with the InPro sensor.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. He will be glad to assist you.

2 Important notes

2.1 Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the InPro 5500 i sensor.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to safety instructions and warnings of potential danger which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies additional information and instructions which, if neglected, could lead to defects, inefficient operation and possible loss of production.

2.2 Intended use

METTLER TOLEDO InPro 5500 i sensors are intended solely for inline measurement of the carbon dioxide concentration in liquids, as described in this instruction manual.

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose.

The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (housings, transmitters, etc).
- observance of all safety regulations governing the equipment in which the sensor is installed.

- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with Mettler-Toledo Process Analytics in the event of any uncertainties.

2.3 Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. housings, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The thermal conductivity sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see spare parts list, "Section 9.3").
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.

2.4 Examples of some typical applications

Below is a list of examples of typical fields of application for the carbon dioxide sensor. This list is not exhaustive.

Measurement in liquids:

- Brewing
- Beverage filtration
- Filling stations for beverage carbonation

3 Product description

3.1 General information

The InPro 5500 i is intended for use in brewery and CSD applications measuring low to high carbon dioxide levels in-line in liquids.

The **carbon dioxide sensor InPro 5500 i** with integrated temperature probe is used for measurement of carbon dioxide from low to high concentration levels (0 ... 15 g/L CO₂).

The sensor is **compatible with CIP/SIP** (Clean-In-Place/Steam-in-place).

The InPro 5500 i sensor with ISM functionality offers Plug and Measure as well as enhanced diagnostics features.

3.2 Principle

The InPro 5500 i utilizes the thermal conductivity principle for determining CO₂ concentration. Thermal conductance is a measure of the quantity of heat that passes in unit time through a plate of particular area and thickness when its opposite face differs in temperature by one kelvin.

3.2.1 Thermal conductivity measurement of CO₂ in liquids

Each gas has a specific thermal conductivity. Therefore, the system is able to measure a gas concentration in binary mixture when this gas has a clearly distinct thermal conductivity to that of the background gas. The binary background gas mixture used to measure CO₂ for the InPro 5500 i sensor is air (N₂/O₂).

The InPro 5500 i uses the technique of gas diffusion through a membrane body that is immersed into the fluid to be measured (beer or carbonated soft drink) The diffused CO₂ gas is then measured by the incorporated thermal conductivity chip.

The sensor head space volume containing the TC-sensor is purged periodically with air (purge gas). The used purge gas builds the base line for the CO₂ measurement and makes the measurement highly CO₂ selective. As CO₂ diffuses through the membrane, the thermal conductivity rises. Using air as the background gas eliminates interference from small amounts of N₂/O₂ being present in the measured process liquid.

3.2.2 Determination of CO₂ concentration

The rise factor (RAMP) the diffusing gas creates on the TC-chip (resulting voltage) is proportional to the partial pressure of the CO₂ in the process liquid $p(\text{CO}_2)$.

Using its integrated temperature sensor, the InPro 5500i directly measures the process temperature for accurate temperature compensation. The temperature compensation is needed due to the dependency of the CO₂ permeability of the gas selective membrane.

The measured partial pressure of CO₂ and the process temperature allow the exact calculation of CO₂ concentration (g/L) using Henry's law.

Henry's Law:

At a constant temperature, the amount of a given gas that dissolves in a given type and volume of liquid is directly proportional to the partial pressure of that gas in equilibrium with that liquid.

$$c(\text{CO}_2) = p(\text{CO}_2) * H$$

$p(\text{CO}_2)$ = partial pressure CO₂

$c(\text{CO}_2)$ = concentration CO₂ (g/L)

H = Henry's constant

3.2.3 Measurement cycles with the InPro 5500 i

Purge cycle (app. 10 sec)

The purge gas flushes the measuring chamber until a flat base line is determined by the TC-Chip.

CO₂ measuring cycle (app. 10 sec)

The diffusing CO₂ increases the measured voltage (thermal conductivity) and the digital sensor calculates the CO₂ concentration.

The calculated value is sent digitally to the connected transmitter.

ISM sensor:

The InPro 5500 i is equipped with Intelligent Sensor Management (ISM).

Principle: The sensor head contains a microprocessor that monitors and controls the sensor and also stores all sensor data. The data is accessed via the transmitter.

The following data are stored permanently in the sensor:

- type of sensor
- serial no.
- software version
- hardware version
- part no.
- operating time
- calibration time and calibration date
- calibration table/data
- CIP/SIP cycles

3.3 Scope of delivery

Each sensor is supplied fully assembled, factory-tested and calibrated for correct function together with:

- a quality control certificate
- factory calibration certificate is stored in the sensor's electronics
- inspection certificates 3.1 (complying with EN 10204)
- protective cap for safe storage

Note: Make sure the sensor is stored in a dry environment when not in use.-

Note: If the MembraCap (membrane cap) is removed or replaced the factory calibration is void. The sensor then needs recalibration.


3.4 Equipment features

The InPro 5500 i is available with three different process connections for CO₂ measurement in brewery and carbonated soft drink applications.

- Varivent DN40
- TriClamp 2"
- 28 mm; M42

The InPro 5500 i can be used in applications where CIP and SIP cycles are applied. The sensor has a protective function where the CO₂-measurement loop shuts off during elevated temperatures (> 55 °C) and turns back on when it cools down (< 50 °C). This function ensures that the TC-chip will not be damaged during CIP/SIP.

Picture:

4 Installation

The digital Plug and Measure feature allows the user to measure CO₂ levels immediately after installation.

The sensor is recognized automatically by the transmitter and all relevant data is sent to the transmitter for automatic configuration. The CO₂ level is then displayed immediately.

The sensor is factory calibrated and the data is stored in the sensor. Calibration data therefore does not need to be entered by the user (to achieve best performance at given temperature we recommend to perform a process calibration).

4.1 Mounting the sensor

Important! Remove the protection cap before mounting the sensor.

Do not remove the MembraCap (membrane cap) that is screwed onto the sensor. If this is done the sensor needs recalibration.

4.1.1 Mounting and positioning the sensor

Attention! Do not turn the sensor counter-clockwise when installed in a 28 mm housing. The Membra-Cap might be loosened.

The sensor is available in three versions and can be fitted without a housing to the following process connections/armatures:

- Sanitary DN40 VARINLINE® Type N Connection
- 2" Tri-Clamp® Sanitary Connector
- 28 mm M42 weld on sockets
- 28 mm M42 VARINLINE® stationary adapters

The sensor must be installed in a way that the MembraCap is in contact with the process fluid to be measured.

- Make sure the sensor is mounted from the side or from the top
- Avoid mounting the sensor from the bottom
- Mount the sensor in a flow ascending position
- Avoid a position near to pump discharges

Flow requirement:

For the sensor's proper functioning the flow in the pipe needs to achieve a min. flow of 0.5 m/s

4.1.2 Removing the sensor and storage

- To remove the sensor from the process the pipe needs to be drained. Avoid counter-clockwise rotation in case where the sensor is installed on a 28 mm socket.

- Dry the MembraCap with a soft tissue
- Stop the purge air flow
- disconnect the purge air connection
- disconnect the 5-pin connector
- Refit the protection cap and store the sensor in a dry environment

4.2 Purge air supply

The sensor needs a purge air supply for correct functioning. The purge gas also keeps the measuring gas chamber dry and prevents humidity build up that could harm the sensor.

Note: If the purge air flow is unstable or interrupted by any means the InPro 5500 i will deliver a fault message to the transmitter and measurement circuit will be shut off for sensor protection.

4.2.1 Requirements of the purge gas supply:

Any air supply that fulfills the requirements below is suitable to be used with this sensor:

- delivered pressure: min. 1 – 2 barg
- particulate free (40 micrometer)
- condensate free (dew point below – 10 °C)

Note: Since the purge air supply also protects the sensor from humidity it needs to be functioning at all times, even if the sensor's measurement is not needed.

If the transmitter shows an alarm the following needs to be checked.

Possible source of unstable or interrupted air supply:

- contaminated filter on the air supply conditioner
- plugging of the outlet tube of the sensor
- insufficient inlet pressure

4.2.2 Connecting the purge air to the sensor

To connect the purge gas supply to the sensor use a flexible nylon, PVC or a stainless steel hose.

The sensor is supplied with a 6mm Swagelock connection for this purpose.

In cases where the air supply requirements cannot be fulfilled METTLER TOLEDO offers a simple purge air conditioner that can be mounted between the air supply and the sensor (Part number: 30 034 319).


4.3 Connections

4.2.1 Connecting the sensor to an RS485 cable

 The sensor is connected to the transmitter via an RS485 5-pin data cable. The cable ensures a secure connection between the transmitter and the sensor under harsh industrial conditions.

 To connect the data cable to the sensor align the slit of the connector with the pins in the plug. Then tightly screw the plug to fasten the two parts.

4.2.2 Connecting the cable to the transmitter

 **Note:** Cable assignment can be found in the METTLER TOLEDO cable instruction manual.

 **Note:** For connecting the cable to the terminals of the transmitter, please also refer to the instructions given in the METTLER TOLEDO transmitter manual.

The cables are available from METTLER TOLEDO in various lengths:

Connect the data cable to the transmitter as described in the tables below.

RS 485 cable		Transmitter		
		M400	M800	
			Channel 2	Channel 4
Color	Function	TB4	TB2	TB4
brown	24 DC+	1	9	9
black	24 DC-	2	10	10
gray	shield	6	12	12
yellow	shield	6	15	15
blue	RS485-	7	13	13
white	RS485+	8	14	14

5 Operation

5.1 Start-up

Each sensor is supplied ready to use. Before using remove the protecting cap.

5.2 Configuration

5.2.1 Sensor detection

Before installing the CO₂ sensor, please refer to the manual for the transmitter and configure the transmitter for automatic sensor detection. The M800 transmitter does not need any pre-configuration.

5.2.2 CO₂ solubility and temp. factor

The sensor is delivered factory calibrated and is set up to measure in beer as the default.

The sensor offers a choice of CO₂ solubility for measurement in beer, cola or water. For other beverages the user has the possibility to enter individual values for CO₂ solubility and temperature factor.

Default values for measurement in beer (Valid for temperatures –5 ... 50 °C):

CO₂ solubility (A): 1.420 g/L

Temp. factor (B): 2485

Values for pure water:

CO₂ solubility (A): 1.471 g/L

Temp. factor (B): 2491

Values for cola:

CO₂ solubility (A): 1.345 g/L

Temp. factor (B): 2370

For beverages where the user knows the exact CO₂ solubility and the temperature factor the values can be changed individually. If the user desires to evaluate the solubility and temperature factors they can be evaluated with the following formulas.

$$H_{CO_2} = A * \exp (B * (1/T - 1/298.15))$$

$$C_{CO_2} = H_{CO_2} * p_{CO_2}$$

H_{CO₂} : Calculated CO₂ solubility (Henry constant) at measured process temp.

A: Solubility of CO₂ (g/L at 25 °C)

B: Temperature factor (valid for –5–50 °C)

C_{CO₂}: Calculated CO₂ concentration

5.3 Calibration

5.3.1 Purpose of calibration

Information about the calibration, you find also in the manual of the transmitter.

Calibration should be performed at least after each change of the MembraCap.

Since the measurement of the thermal conductivity is dependent on the diffusional properties of the used MembraCap, a calibration after exchange is mandatory to regain accurate CO₂ measurements.

5.3.3 Single point calibration in-line, off-line

The sensor is designed to provide accurate CO₂ measurement with a simple slope calibration using a pure CO₂ gas or using a gas mixture where the partial pressure of the CO₂ content is known.

The calibration can be performed in-line or by dismounting the sensor and performing the calibration in the workshop.

METTLER TOLEDO offers a calbox (P/N: 52 300 400) as a support for the sensor for workshop calibration. The CalBox has a 6 mm hose connection to connect the calibration gas.

Note: If the calibration is done in the workshop, make sure a purge air source is used.

Note: The air supply (purge air) needs to be active at all times while performing a calibration.

Calibration procedure:

- Start the calibration by pressing the cal button on the transmitter and choosing "1-point calibration"
- Process pipe needs to be filled with the calibration gas
- After exposing the sensor to the calibration gas confirm on the transmitter that the "Sensor is in Gas"
- Enter the known CO₂ partial pressure in the calibration menu in mbar or hPa and continue
- Per default the sensor starts calibration when CO₂ measurement is stabilized (Auto-CAL Mode)
- For user convenience the Auto-Cal mode can be changed to manual cal mode. In this case the calibration can be started as soon as the sensor delivers a measured value (approx. every 20 sec.), provided a stable reading is given
- Once the calibration is done, save the calibration onto the sensor by selecting "Save Adjust"
- A message "calibration successful" appears and the user is asked to "Re-install sensor"
- Now the sensor is recalibrated and ready for use

-  **Note:** If the message "Cal out of Limit" appears please check the gas supplies (calibration gas and purge gas) and/or if the MembraCap needs to be exchanged. In this case please abort the calibration.
-  **Note:** If only a verification of the calibration is needed, choose "safe calibrate" in the transmitter menu. There is also the choice for aborting the calibration "save abort".
-  **Note:** When the MembraCap is exchanged the sensor needs to be operated at least 15 min. before performing a calibration.
-  **Note:** Be aware that when calibration is done with CO₂ reference gas the sensor will show slightly different CO₂ values in gas. The reason for this is that the sensor's electronic use a compensation to give accurate measurement in liquids. For calibration verification purposes with reference gas please follow the steps in chapter 6, comparing the new and old slope values.

5.3.4 In-line process single-point calibration

If the user requires that the sensor measures the same values as with established off-line reference CO₂ measuring methods, the InPro 5500 i offers a single-point process calibration option.

Calibration procedure:

- Start calibration by pressing the cal button on the transmitter and select CO₂ calibration
Type = Process
- Take a sample from the process fluid at the same time from which the CO₂ value has to be determined by an established reference method and confirm capture
- After analyzing the CO₂ value go back to the transmitter and start cal again and enter the reference value in mg/L, V/V, mbar or hPa
- Once the calibration is done, save the calibration to the sensor by confirming "Save Adjust"
- A message "calibration successful" appears
- Now the sensor is process calibrated and ready for use

-  **Note:** If message "Cal out of Limit" appears please check your reference measurement or if the membrane cap needs to be exchanged. In this case please abort the calibration.

6 Maintenance

6.1 Inspection of the sensor

6.1.1 Visual inspection

To check your sensor, we recommend the following procedure:

- The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Before calibration always examine the Membra-Cap optically for signs of damage. The Membra-Cap must be intact and clean. Dirty membranes should be wiped clean using a soft, moist tissue.
- The membrane body has to be replaced once to twice per year to maintain max. measuring performance.
- The MembraCap needs also to be replaced if the sensor cannot be calibrated or the MembraCap shows sign of mechanical damage.
- The Intelligent Sensor Management (ISM) Dynamic Lifetime Indicator will tell the user when a MembraCap change needs to be performed.

Attention! Do not use any aggressive cleaning agents. This could damage the sensor or lead to faulty measurement. If the dirt can not be removed with a moist tissue, use gently non-aggressive solvent (alcohol). Any other chemicals shall not be used. This could damage the sensor or lead to a faulty measurement.

6.1.3 Testing the sensor with a transmitter

For verifying the sensor, check measurements vs. reference (or known) sample values.

This can be done in-line in process solution or using CO₂ gas where the partial pressure is known. For off-line/workshop verification only use a calibration gas of a known value.

If the measured values differ from the desired values, a new process or a calibration with reference gas might be performed.

Note: When the MembraCap is exchanged the sensor needs to be operated at least 15 min. before performing a calibration.

If the sensor does not perform as expected after replacing the MembraCap send the sensor to your local METTLER TOLEDO representative for inspection.

Example:
Sensor verification using a transmitter and a calibration gas:

- Press Cal button “Calibrate Sensor” and confirm.
- Chose “1-point CO₂” calibration and confirm.
- Start calibration when sensor is in gas.
- Enter the partial pressure of the calibration gas.
- The transmitter will now show the slope $S = **mV$ and the baseline values $BL = **mV$ and asks for “Save adjust”.
- If you want the sensor to measure based on the new calibration data choose “Save adjust”.
- If you want only to save the calibration data chosse “Save calibrate” (the sensor in this case still measures with the previous adjusted calibration data).
- for verification purposes with a calibration gas the slope values can be compared with the previous calibration values. The previous 3 calibrations are saved in the ISM-Info menu (Info Menu: Cal Data)

Note: The sensor now still uses the previous calibration data for the measurements. The sensor will only measure with the new calibration data by using the command “save adjust”. “Safe calibrate” stores the calibration data for verification tracking purposes only.

Note: If the notice “cal out of limit” appears, please check the purge gas supply and calibration gas supply. If you still get the notiice, replace the Membra-Cap and perform a new calibration.

6.1.4 ISM diagnostics

Dynamic Lifetime Indicator: DLI

The DLI provides information about the remaining lifetime of the MembraCap. As long as the DLI is above zero days the system is within the specified accuracy after a calibration. If the DLI is zero the MembraCap shall be replaced. Contributing factors for wear of the MembraCap are:

- temperature during measurement
- number of CIP cycles
- number of SIP cycles

DLI calculation: Using the above parameters actual sensor stress is calculated. The accumulated sensor load divided by the elapsed time is the basis of the calculation of the remaining lifetime.

Adaptive Calibration Timer: ACT

The ACT provides information as to when the next calibration is required to ensure measurements will remain within the specified accuracy. This calculation is based on the sensor wear information.

Calibration history

The last three calibrations and the factory calibration data are stored in the sensor memory. These data can be read out with a transmitter using the "Info Menu".

The calibration history gives valuable information regarding the quality of the calibration and the wear of the MembraCap.

6.1.5 ISM design

The integrated ISM functionality allows extensive monitoring of the sensor. The following parameters are stored in the sensor:

- serial no.
- type of sensor
- calibration data
- CIP/SIP counter
- slope
- baseline

Whenever the sensor is connected to a transmitter the following automatic test procedures are implemented:

- digital communication
- Plug and Measure
- pre-calibration
- predictive maintenance

6.2 Changing the MembraCap

The sensor has been designed to allow the user to exchange the service parts within the shortest time frame.

This is achieved by reducing the maintenance to only one part, the replacement of the MembraCap.

The MembraCap consists of one piece and incorporates a polymeric membrane that controls the diffusion of CO₂ to the measuring chamber where the thermal conductivity is measured.

Replacement procedure:

Attention! When the MembraCap is detached, take care not to damage the inner parts of the sensor. Damage and soiling of the inner parts, especially the thermal conductivity chip, may influence the signal or destroy the sensor.

When changing the MembraCap, please observe the following instructions:

Attention! Make sure that this maintenance step is carried out in a clean environment.

Attention! Make sure you hold the sensor upside down (Fig. 1, 2 and 3). Otherwise the diffusor might fall out and the thermal conductivity chip will be exposed. If the diffusor falls out clean it gently with alcohol and place it back on top of the sensor.

Fig. 1: Hold sensor in upside down position during maintenance e.g. MembraCap change

Fig. 2: Sensor view with removed diffusor

Fig. 3: Sensor with removed MembraCap

1. Dry sensor carefully with a soft tissue
2. Unscrew the MembraCap sleeve from the sensor shaft and carefully pull it off the sensor (Fig. 4).
3. Place the new MembraCap on top of the inner body of the sensor shaft (Fig. 5).
4. Carefully slip the MembraCap over the sensor and screw it down as tightly as possible by hand (do not use a tool) [Fig. 6 to 7].
5. After each exchange of the MembraCap, the DLI must be reset manually using the transmitter.
6. After each exchange of the MembraCap, the sensor has to be recalibrated by a 1-point calibration using a CO₂ calibration gas or a process calibration in measuring solution. Make sure to reset DLI before performing a calibration. The sensor will automatically prompt for calibration.

Fig 4: Unscrewing MembraCap counterclockwise

Fig. 5: Placing the new MembraCap on top of the sensor

Fig. 6: Screwing on the new MembraCap

Fig. 7: Sensor with newly mounted MembraCap, ready for calibration

Attention! The quality of this calibration is critical for the sensor performance and the accuracy of the diagnostics.

6.3 Error Handling

Message displayed in SERVICE menu	Message in INFO menu	Potential cause	Recommendations
ChX error CO ₂ slope > errorlimit high	ChX error CO ₂ slope > errorlimit high	<ul style="list-style-type: none"> aged membrane faulty cal data input 	<ul style="list-style-type: none"> change membrane enter correct cal data
ChX error CO ₂ slope < errorlimit low	ChX error CO ₂ slope < errorlimit low	<ul style="list-style-type: none"> aged membrane low purge gas flow faulty cal data input 	<ul style="list-style-type: none"> change membrane assure correct purge gas flow enter correct cal data
ChX CO ₂ overrange	ChX CO ₂ out of range ChX CO ₂ not reliable	<ul style="list-style-type: none"> leak membrane insufficient purge gas flow humidity in sensor 	<ul style="list-style-type: none"> change membrane assure correct purge gas flow contact Mettler-Toledo service
ChX CO ₂ underrange	ChX CO ₂ out of range ChX CO ₂ not reliable	<ul style="list-style-type: none"> humidity in sensor drift in electronics 	<ul style="list-style-type: none"> contact Mettler-Toledo service contact Mettler-Toledo service
ChX Temp. overrange	ChX Temp. out of range ChX CO ₂ not reliable	<ul style="list-style-type: none"> temp. above measuring range (e.g. during CIP/SIP) 	<ul style="list-style-type: none"> decrease process temperature no further action required
ChX Temp. underrange	ChX Temp. out of range ChX CO ₂ not reliable	<ul style="list-style-type: none"> temp. below measuring range 	<ul style="list-style-type: none"> increase process temperature
ChX CIP counter expired	ChX CIP counter expired	<ul style="list-style-type: none"> too many CIP's 	<ul style="list-style-type: none"> NA
ChX SIP counter expired	ChX SIP counter expired	<ul style="list-style-type: none"> too many SIP's 	<ul style="list-style-type: none"> NA
ChX TC sensor cut off	ChX CO ₂ not reliable	<ul style="list-style-type: none"> temp. above measuring range TC sensor failed humidity in sensor 	<ul style="list-style-type: none"> decrease process temperature check sensor for humidity or membrane defects ▶ reset CO₂ measurement in transmitter menu
ChX board error	ChX hardware error	<ul style="list-style-type: none"> defectuous el. board 	<ul style="list-style-type: none"> contact Mettler-Toledo service
CHX TC-sensor fail	ChX hardware error	<ul style="list-style-type: none"> leak membrane humidity in sensor 	<ul style="list-style-type: none"> check sensor for humidity or membrane defects ▶ reset CO₂ measurement in transmitter menu
CHX valve open	ChX CO ₂ not reliable	comes up if: <ul style="list-style-type: none"> TC sensor fails membrane error active temp. above measuring range 	<ul style="list-style-type: none"> NA
ChX SW error	ChX software error	<ul style="list-style-type: none"> FW of measuring and communication processor not compatible 	<ul style="list-style-type: none"> contact Mettler-Toledo service
ChX error CO ₂ slope > warninglimit high	ChX error CO ₂ slope > warninglimit high	<ul style="list-style-type: none"> aged membrane faulty cal data input 	<ul style="list-style-type: none"> change membrane enter correct cal data
ChX error CO ₂ slope < warninglimit low	ChX error CO ₂ slope < warninglimit low	<ul style="list-style-type: none"> aged membrane low purge gas flow faulty cal data input 	<ul style="list-style-type: none"> change membrane assure correct purge gas flow enter correct cal data

ChX BL out of range	ChX hardware error	<ul style="list-style-type: none"> insufficient purge gas flow humidity in sensor drift in electronics 	<ul style="list-style-type: none"> assure correct purge gas flow contact Mettler-Toledo service contact Mettler-Toledo service
ChX membrane error	ChX change membrane	<ul style="list-style-type: none"> leak/defectuous membrane 	<ul style="list-style-type: none"> change membrane
ChX DLI expired	ChX DLI expired	<ul style="list-style-type: none"> membrane life time potentially expired 	<ul style="list-style-type: none"> change membrane reset DLI counter
ChX Cal required	ChX Cal required	<ul style="list-style-type: none"> potentially too long time window since last calibration 	<ul style="list-style-type: none"> execute calibration
ChX Temp. not reliable	ChX CO ₂ not reliable	<ul style="list-style-type: none"> short or interrupt in temp. line defectuous temp. sensor defect on el. board 	<ul style="list-style-type: none"> contact Mettler-Toledo service contact Mettler-Toledo service contact Mettler-Toledo service

7 Storage

For storage the sensor should be clean and dry. The protection caps have to be placed on the sensor and the cable connectors as soon as the sensor is not in service.

The protection cap of the sensor incorporates a desiccant capsule to keep the sensor dry.

Make sure to use a fresh capsule before mounting the protection cap to the sensor

Attention! The sensor must not be stored without the protection caps.

8 Product specification

8.1 Certificates

Each sensor is delivered with a set of **3.1 certificates** (complying with EN 10204).

All wetted metal parts (sensor shaft, MembraCap are identified with an engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part (sensor shaft and membrane body) is polished in order to achieve a surface roughness lower than 0.4 µm (16 µin). This represents a roughness grade number of N5 (according to ISO 1320:1992).

8.2 Specifications

InPro 5500 i	
Measurement principle	Thermal conductivity
Working conditions	
Permissible pressure range during measurement	0 ... 20 bar absolute (except pressure spikes)
Mechanical pressure resistance	= permissible pressure range
Permissible temperature range during measurement	-5 ... 50 °C [23 ... 122 °F]
Mechanical temperature range	-5 ... 120 °C [23 ... 248 °F]
Sensor performance	
Operation range	0 ... 10 bar p(CO ₂) 0 ... 15 g/L CO ₂ 0 ... 7 V/V CO ₂
Accuracy in fluids	± 1 % within ± 5 °C of calibration temperature ± 2 % over temperature range 0 ... 50 °C
Measurement cycle	< 20 s
Flow requirements	min. 0.5 m/s
Design features	
Process connections	– Varivent™ Type N – Tri-Clamp® 2" – 28 mm with cap nut M42
Temperature compensation	automatic with built-in RTD
Cable connection (digital)	5 pin, RS-485 data cable
O-ring material (wetted parts)	EPDM (other material on request)
CO ₂ selective membrane material	PTFE/silicone (reinforced with steel mesh)
Material sensor body (wetted parts)	316L stainless steel
Surface roughness of wetted metal parts	N5 (Ra = 0.4 µm [16 µin])
IP protection class	IP 67
Certification MaxCert™	
Material certificate 3.1	Yes
Final inspection certificate	Yes

9 Ordering information

For more detailed information refer to the technical data sheet. Ask your local distributor.

9.1 Sensors

Sensors	Order No.
InPro 5500i/Varivent Type N	30 034 265
InPro 5500i/Tri-Clamp 2"	30 034 266
InPro 5500i/28mm/M42	30 034 264

9.2 Accessories

Accessories	Order No.
CalBox	52 300 400
Purge gas conditioner	30 034 319

Cables

Data cable (5 pin) for InPro 5500i
temperature range -30 ... 80 °C (-22 ... 176 °F)

RS485 data cable 2 m	52 300 379
RS485 data cable 5 m	52 300 380
RS485 data cable 10 m	52 300 381
RS485 data cable 15 m	52 206 422
RS485 data cable 25 m	52 206 529
RS485 data cable 50 m	52 206 530

9.3 Spare parts

Spare parts	Order No.
MembraCap	30 034 318
O-ring set (Varivent)	30 038 616

9.4 Recommended transmitters

Transmitters	Order No.
M400, type 3	52 121 350
M800, 2-channel	52 121 813

Note: No housings is needed, the sensor is designed to fit directly into the three offered process connections!

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- CH** **Mettler-Toledo (Schweiz) GmbH**, Im Langacher, Postfach, CH-8606 Greifensee
- D** **Mettler-Toledo GmbH**, Prozeßanalytik, Ockerweg 3, D-35396 Gießen
- F** **Mettler-Toledo, Analyse Industrielle S.A.S.**, 30, Boulevard Douaumont, F-75017 Paris