
METTLER TOLEDO AT-SE Kit

Mounting instructions
Conversion kit 210680
starting page 3

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METTLER TOLEDO

What AT-SE represents:

The AT-SE Kit can be used to convert every METTLER TOLEDO AT balance to allow the evaluation electronics to be set up separately (SE) from the weighing unit (max. distance = 5 m).

What is this conversion used for?

The conversion allows the weighing cell to be used at exposed sites with the evaluation unit (abbreviated to EU) located outside this zone and thus protected. Data are transmitted by means of the connection cable.

The conversion kit contains all parts needed to convert an AT balance into an AT-SE balance.

The parts list at the end includes all parts needed for the conversion.

To ensure that no fastening components (screws etc.) are missing, all such components are included with the conversion kit although, for instance, the bolts used to fasten the balance board could be reused for installation in the EU after board removal.

The conversion must be performed by a METTLER TOLEDO dealer.

Conversion principle:

The evaluation electronics of an AT balance are disassembled and built into a new, separate housing (EU).

The sequence of the work steps is listed in the following description "Converting the balance". The item numbers in brackets mentioned in these conversion instructions match the item numbers of the spare parts list/exploded view diagram of the AT service manual, 704008.

The parts list on page 6 specifies the installation position of each part and serves as an additional aid.

Spare parts for the AT-SE

The spare parts for the AT-SE are listed in the AT service manual.

Converting the AT balance to the AT-SE

- Mount the three rubber feet (screw feet) with the countersunk screws M4x10 and the nuts (for locking) on the housing base of the EU.
- Secure the connector bracket to the 20-pin EU cable connector (display) using screws M2.5x6.
- Fasten the connector bracket with screws M3x5 from above to the housing base.
- Fasten the second, 25-pin EU cable connector with the assembly bolts, washers and nuts to the housing base at the back. (Connector inside the housing base, attach assembly bolts from outside.)
- Fasten display console to the housing base at the front with screws M4x8.
- Remove display unit from the AT balance and fasten to EU display console with screw M4x8.
- Remove rear housing panel (33) of AT balance.
- Pull off housing cover (15) together with the upper sliding door (27) backward. Ensure that the guide rails (63) + (64) do not project from the guide otherwise the rollers may drop out.
- Remove interface board (42) after disconnecting the cable connection at the CPU board.
- Remove CPU board (37) after undoing the cable connections.
- The CPU board must be fastened to the EU housing base using the board supports in the conversion kit *1).
- Plug EU cable connections into CPU board.
- Mount interface board in the EU housing base on the board supports and fasten with screws M3x5 *2).
- Plug connector of the interface into CPU board.
- Mount EU cover (top housing) on the EU housing base (the Mettler logo should be at the rear on the right) and fasten from below with screws M4x8.
- Insert cassette cover at EU rear panel.
- Install weighing cell board in AT balance and fasten with the 4 special screws (with collar). The board must be capable of being moved to and fro.

- Plug the 3 cable connections into cell board.
- Insert housing cover together with the upper sliding door.
- Fasten the AT-SE rear panel with screws M4x12.
- Screw the 25-pin connector of the cell board using the assembly bolts to the AT-SE rear panel. Nuts and washers are not needed!
- Remove AT front frame (front glass) (72).
- Fasten AT-SE front frame and cover with screws M4x16 to the balance.
- Attach connection cable (length = 5 m) to the balance and EU.
- Connect AC adapter.
- Check functioning of the AT-SE.

*1) It must be ensured that the connector shielding has good mechanical contact with the rear panel of the housing base. If not, bend the contact spring appropriately.

*2) It must be ensured that the connector shieldings have good mechanical contact the rear panel of the housing base. If not, bend the contact springs appropriately.

Parts list/standard equipment for the AT-SE conversion

Part	Designation	Installation site
1	Housing base, cpl.	EU (evaluation unit) base
1	Weigh. cell board AT-SE	AT-SE balance
1	Display console	EU housing base, front
1	Cover	AT-SE front frame (display cover)
1	EU cable	20-pin connector → connector bracket, 25-pin con. → housing base rear panel
2	Cheese H. screw M2.5x6	20-pin connector at connector bracket
6	Cheese H. screw M3x5	Connector bracket → housing base 2 / interface board → board supports 4
7	Cheese H. screw M4x8	Display console 2 / display 1 / cover on EU housing base 4
3	Cheese H. screw M4x12	AT-SE rear panel
2	Cheese H. screw M4x16	Cover → AT-SE front frame
3	P. Countersunk scr. M4x10	Rubber feet → housing base from below
3	6KT-05MU M4	Rubber feet → housing base from above (inside)
3	Rubber feet	Housing base from below (screw feet)
2	Bolt set	25-pin connector (EU cable) → housing base rear panel 25-pin connector (cell board AT-SE) → balance rear panel
1	Connector bracket	Housing base
4	Board supports	CPU board on housing base
4	Special screws	Weighing cell board AT-SE
1	Cassette cover	EU rear panel
1	Cover, cpl.	EU housing cover (screwed to housing base)
1	Front frame SE, cpl.	AT-SE balance
1	Rear panel	AT-SE balance
1	Connection cable (5 m)	AT-SE to EU
1	Mounting instructions	

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© Mettler-Toledo GmbH 1999 704267 Printed in Switzerland 9908/33.31
Mettler-Toledo GmbH, Laboratory & Weighing Technologies, CH-8606 Greifensee, Switzerland