
**Operating Instructions
Bedienungsanleitung
Mode d'emploi**

**METTLER TOLEDO
GA42 Printer
Drucker GA42
Imprimante GA42**

METTLER TOLEDO

Overview, Übersicht, Sommaire

Introduction	3
Safety measures	3
Startup procedure	4
Operation	5
Printer test	5
Maintenance	6
Malfunctions	6
Cleaning	7
Accessories, Spare parts	7
Technical data	8
CE declarations of conformity	24
Safety tests	25
ISO 9001 certificate	25
Explanation of the validation	26
Product development process	27

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Introduction

The GA42 is a rugged and simple to operate dot matrix printer for METTLER TOLEDO titrators, pH meters, density meters, refractometers, and melting point apparatus, as well as for METTLER TOLEDO balances with an RS232C interface (balance lines: AB, B, PB, CB, GB, PG-S). The printout on normal paper fulfills the requirements of modern quality assurance systems (GLP, GMP, ISO 9000, etc.). These define, among other things, the principles and duties for the user of measuring instruments or inspection, measuring and test equipment. The GA42 complies with the principle of traceability of measurements by precisely documenting the origin of data, as well as calibration procedures. It fulfills the obligation to archive measurement results by its chemically resistant, thermally stable and lightfast print. On connection to METTLER TOLEDO balances, the GA42 prints all data sent by the balance without restriction, that is not only weight values and calibration procedures, but also the results of all applications, such as piece counting, built into the balance.

Safety measures

The GA42 has been tested for the connections and intended purposes documented in these Operating Instructions. However, this does not absolve you from the responsibility of performing your own tests of the product supplied by us regarding its suitability for the methods and purposes you intend to use it for. You should therefore observe the following safety measures.

Risk of electric shock

- Never open the housing of the printer yourself! Have the printer serviced only by METTLER TOLEDO service!

Risk of explosion


- Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight (explosion hazard due to spark formation, corrosion caused by the ingress of gases).

Caution


- Wipe off splashed liquid immediately! The printer housing is not watertight.
- Follow the cleaning instructions on page 7!

Startup procedure

Configuring the interface

- Open cover of paper compartment (1) and remove paper roll.
- Press plastic tab (2) backward and open the battery compartment cover.
- Leave DIP switch 1 in the OFF position! It has no meaning.
- Set DIP switch 2 of the printer accordingly.
 - ON: For DL36 KF Coulometer, DL31 / 38, DL50 Graphix, DL50 / 53 / 55 / 58 Titrators
MP125, MA130, MC126, MO128 portable Meters
DELTA 340 / 345 / 350 / 355 pH Meters
DA-100M / 300M / 310M, DE40 / 50 / 51 Density Meters
RE20B / 40 / 50 Refractometers
 - OFF: For DL18 / 20 / 21 / 25 / 35 / 40 / 40RC / 40GP Titrators
QUATTRO MP225, MP230, MA235, MC226, MPC227
RA-510M Refractometer
FP62 Melting Point Instrument and FP800 / 800HT Thermosystems
all METTLER TOLEDO balances with an RS232C interface


Inserting paper

- Feed paper through the slot in the printing unit in the direction of the arrow.
- Switch on printer with toggle switch at the rear.
- Press and hold  until sufficient paper has been fed in.

Attaching the printer

- Attach connection cable of the GA42 to the printer interface of the instrument and to the printer. After configuring the attached instrument for the GA42, the printer is ready for operation. Order numbers for connection cables, see page 7.

Operation

- Pressing the printer key  advances the paper 1 line (line feed).

On connection of the GA42 to METTLER TOLEDO balances:

- Pressing the balance key with the  symbol prints out the stable weight value.
As long as the balance does not show a stable weight value, nothing will be printed out.
With PR- and SR-balances, an automatic printout is possible. The printout includes weight with time, date and sample number.

Printer test

The GA42 is equipped with a self-test which automatically prints out the character set of the printer.


Starting the test

- Switch on GA42 and at the same time press and hold  for a few seconds.

Terminating the test

- Switch off the GA42.

Maintenance


Changing paper

- Insertion of a new roll of paper is described on page 4.
- Order number for paper, see page 7.

Replacing the ribbon

- Pull paper out of printing unit (1).
- Remove ribbon cover (2, 3, 4).
- To remove ribbon cartridge, press in direction of arrow (5).
- Insert new ribbon and if need be, take up slack with tensioning wheel (6).
- Refit ribbon cover.
- Insert paper and feed through slit in the printing unit as far as it will go.
- Press and hold  until sufficient paper has been fed in.

Malfunctions

Malfunction / Message

Printout unreadable

Green pilot lamp does not light up
Green pilot lamp flashes slowly


Possible cause

- Ribbon worn / entangled
- Service life of printing unit at an end
- No power, fuse defective
- Balance waiting for stability / no communication with attached instrument

Rectification

- Change / tension ribbon
- Please contact METTLER TOLEDO service
- Please contact METTLER TOLEDO service

Cleaning

As the printer housing is made of top grade, resistant materials, all commercially available cleaning agents may be used.

Accessories, Spare parts, Consumables

	Type / Order number	Standard
Paper roll, set of 5	72456	2
Ribbon cartridge, black, minimum order quantity: 2	65975	1
Connection cable for DL18 / 20 / 21 / 25 / 35 / 40 / 40RC / 40GP, FP62	17842	–
Connection cable for RA-510M, DA-100M / 300M / 310M	51328094	–
Connection cable for FP800 / 800HT	17859	–
Connection cable for DL50 Graphix, DL50 / 53 / 55 / 58	201508	–
Connection cable for DL36, DL31 / 38, DE40 / 50 / 51, RE20B / 40 / 50	51190362	–
Connection cable for QUATTRO MP225, MP230, MA235, MC226, MPC 227	51302125	–
Connection cable for Portables MP125, MA130, MC126, MO128	51302031	–
Connection cable for DELTA 340 / 345 / 350 / 355	51300162	–
Connection cable for METTLER TOLEDO balances with an RS232C interface	11101051	1
Spare printing unit (may be changed only by METTLER TOLEDO service!)	600133	–
Power line fuse 230 V (T100L250V) or 115 V (T200L250V) – (may be changed only by METTLER TOLEDO service!)	commercially available	–

Technical data

Printer

Printing unit	Dot matrix printer 5x9 dots, 24-character line length, Epson IBM character table No. 4
Printing speed	Min. 1 line per second
Ribbon cartridge	Exchangeable, black
Paper roll	Normal paper, 58 x dia. 51 mm, integrated in housing, commercial size

Printer functions

Printing, lightfast and thermally stable (GLP, GMP, ISO 9000), line feed, self-test


Interface

RS232C

DIP switch 2 in OFF position: 2400 baud, even parity, 7 bits, 2 stop bits

DIP switch 2 in ON position: 1200 baud, no parity, 8 bits, 2 stop bits

Connections: D-sub, 9 pin, male

Designed test conditions

Power line voltage, frequency	115 V or 230 V, admissible fluctuations –20% +15%, 50/60 Hz
Power consumption	140 mA or 70 mA
Altitude	4000 m
Temperature range	0...40 °C
Humidity	Max. relative humidity 80% for temperatures to max. 31 °C, linear decreasing to 50% at 40 °C
Installation category	II
Pollution degree	2
Dimensions / Weight	W x D x H = 157 x 210 x 85 mm / net 1.5 kg (incl. paper roll)

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to both Part 15 of the FCC Rules and the radio interference regulations of the Canadian Department of Communications. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

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Declaration of conformity

We, **Mettler-Toledo GmbH**, Im Langacher, CH-8606 Greifensee, Switzerland, declare under our sole responsibility that the product listed below – marked with “CE” – is in conformity with the directives and standards mentioned.

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Greifensee, 1998

Philip Bodmer
General Manager
Business Unit Laboratory

Stefan Ziegler
Marketing Manager
Business Unit Laboratory

Safety testing of our products

The product listed below has passed the required safety tests for electrical products for the approvals and marks mentioned.

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ISO 9001 certificate for METTLER TOLEDO

The Mettler-Toledo GmbH company, Greifensee, was examined and evaluated in 1991 by the Swiss Association for Quality and Management Systems (SQS), and was awarded the ISO 9001 certificate. This certifies that Mettler-Toledo GmbH, Greifensee, has a quality management system that conforms with the international standards of the ISO 9000 series.

Repeat audits are carried out by the SQS at intervals to check that the quality management system is operated in the proper manner and is continuously updated in relation to changes brought about.

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Explanation of the validation of lab balances and their optional equipment

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All products of the business unit Lab from METTLER TOLEDO including internal software and optional equipment are developed and manufactured within the framework of the quality management system certified in compliance with ISO 9001. They include the following models:

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Precision balances BB, PB, PG, PG-S, PJ, PM, PR

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HG, HR, LJ, LP

Design validation is according to in-house standards for the product and software development process (see schematic representation).

The appropriate validation reports of these process procedures and accompanying documentation such as source codes can be examined after prior agreement with METTLER TOLEDO.

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Philip Bodmer
General Manager
Business Unit Laboratory

Stefan Ziegler
General Manager
Business Unit Laboratory

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To protect your METTLER TOLEDO product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of all METTLER TOLEDO products for years to come. Please send for full details about our attractive terms of service. Thank you.

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© Mettler-Toledo GmbH 1997, 1998, 1999