
User manual



METTLER TOLEDO Compact scales BBA422 / BBA425 / BBK422





Congratulations on choosing the quality and precision of METTLER TOLEDO. Proper use according to this Operating Manual and regular calibration and maintenance by our factory-trained service team ensures dependable and accurate operation, protecting your investment. Contact us about a ServiceXXL agreement tailored to your needs and budget.

We invite you to register your product at www.mt.com/productregistration so we can contact you about enhancements, updates and important notifications concerning your product.

BBA422 / BBA425 / BBK422 Table of contents

Table of contents

		Page
1 1.1 1.2 1.3 1.4	Introduction Safety instructions Description Putting into operation Disposal	5 9
2.1 2.2 2.3 2.4 2.5 2.6 2.7	Operation Switching on and off Zeroing / Zero point correction Simple weighing Weighing with tare Dynamic weighing Printing results Cleaning	11111213
3.1 3.2 3.3 3.4 3.5 3.6 3.7	Settings in the menu Operating the menu Overview Scale settings (SCALE) Application settings (APPLICATION) Terminal settings (TERMINAL) Configuring interfaces (COMMUNICATION) Diagnosis and printing out of the menu settings (DIAGNOS)	15 20 22 23
4 4.1 4.2	Interface description	28
5	Event and error messages	32
6 6.1 6.2	Technical data and accessories Technical data Accessories	34
7 7.1 7.2 7.3	Appendix Information for certified scales in EC countries Safety checks Table of Geo Values	40 40
8	Index	44

Introduction BBA422 / BBA425 / BBK422

1 Introduction

1.1 Safety instructions



CAUTION!

Do not use BBA422 / BBA425 / BBK422 in hazardous areas! Our product range includes special devices for hazardous areas.



CAUTION!

Scales with protection level IP65 are dust-tight and hose-proof to EN 60529. They are suitable for use in dusty environment and brief contact with liquids. Ensure that the scale is dried off again after coming into contact with liquid.

Even with degree of protection IP65 the scale should not be used in environments in which there is a risk of corrosion.

▲ Do not flood the scale or submerge it in liquid.



DANGER!

Electric shock hazard!

▲ Always pull out the mains plug before any work on the device.



DANGER!

Electric shock hazard if the mains cable is damaged!

- ▲ Check the mains cable for damage regularly and replace it immediately if it is damaged.
- ▲ On the rear side of the device, maintain a clearance of at least 3 cm in order to prevent the mains cable bending too much.



CAUTION!

On no account open the device!

The warranty is void if this stipulation is ignored. The device may only be opened by authorized persons.

▲ Call METTLER TOLEDO Service.

BBA422 / BBA425 / BBK422



CAUTION!

Handle the compact scale with care.

The scale is a precision instrument.

- ▲ When the weighing pan has been removed, never clean the area under the load plate holder with a solid object!
- ▲ Do not put excessive loads on the scale.
- Avoid banging the weighing pan.

Note Use with foodstuffs

Parts coming into contact with foodstuffs have smooth surfaces and are easy to clean. The materials used do not splinter and are free of harmful substances.

With foodstuffs, it is recommended to use the supplied protective cover.

- → Clean the protective cover regularly and carefully.
- → Replace damaged or very dirty protective cover immediately.

1.2 Description

This user manual applies to the following types of compact scales:

- Compact scale BBA422... with strain gauge weighing cell, Protection Class IP43
- Compact scale BBA425... with strain gauge weighing cell, Protection Class IP65
- Compact scale BBK422... with MonoBloc, Protection Class IP43

The compact scales are available in a small and large size in various capacities and resolutions.

The power supply is carried out via a built-in power supply device, an internal rechargeable battery with an external mains adapter or an external battery.


One of the following options can also be ordered:

- Additional interface RS232 or RS485
- Ethernet interface
- USB interface
- Digital I/O


Introduction BBA422 / BBA425 / BBK422

1.2.1 Overview

- 1 Display
- 2 Scale specifications
- 3 Load plate
- 4 Adjustable feet
- 5 Keys




- 1 Power supply connection
- 2 Fast and fine pressure equalization, only with Protection Class IP65
- **3** Optional interface
- 4 (Standard) RS interface



BBA422 / BBA425 / BBK422

1.2.2 Display



- 1 Active interface
- 2 Weighing range display
- **3** Battery charge level; only present on scales with a battery
- 4 Symbol for displaying net values
- 5 Symbol for dynamic weighing
- 6 Weight units
- 7 -segment display, 7 digits, with decimal point
- **8** Stability monitor (goes out when a stable weight value is reached)
- 9 Sign
- 10 Identification for changed or calculated weight values, e.g. higher resolution, minimum weight not reached

Introduction BBA422 / BBA425 / BBK422

1.2.3 Keypad

Main functions

Key	Function in operating mode	Function in the menu
0	Switching device on / off, abort	To the last menu item -End-
→0←	Setting scale to zero	Scrolling back
→T←	Taring scale	Scrolling forward
	Transfer key Long key press: Calling up menu	Activating menu item Accepting selected setting

Additional functions

Key	Function
Units	Switching weight unit
Clear	Clear key

BBA422 / BBA425 / BBK422

1.3 Putting into operation


1.3.1 Selecting or changing the location

The correct location is crucial to the accuracy of the weighing results!

→ Select a stable, vibration-free and if possible a horizontal location.










- No direct sunlight
 - No strong drafts
 - No excessive temperature fluctuations







Aligning the scale

Only scales that have been aligned precisely horizontally provide accurate weighing results. The certified scales have a spirit level to simplify alignment.

→ Turn the adjustable feet of the scale until the spirit level's air bubble is inside the inner circle.

Major geographical location changes

The manufacturer adjusts each scale to the local gravity conditions (GEO value). In the event of major geographical location changes, this setting must be adjusted by a service technician. Certified scales must also be recertified observing the national certification regulations. These steps are not necessary for scales with an internal calibration weight.

1.3.2 Connecting the power supply



CAUTION!

Before connecting the scale to the mains, check whether the voltage value printed on the rating plate corresponds with the local mains voltage.

- ▲ Never connect the device if the voltage value printed on the rating plate is different to the local mains voltage.
- → Plug the mains plug into the socket.


 After connection, the device performs a self-test. When the zero display appears, the device is ready to weigh.
- → Calibrate the device in order to obtain the greatest possible precision, see Section 3.3.1.

Introduction BBA422 / BBA425 / BBK422

Note

Partially certified scales (scales with first-level certification) must be certified by an authorized body or by the METTLER TOLEDO Service.

→ Call METTLER TOLEDO Service.



Scales with a built-in battery can work independently from the mains for approximately 30 hours in normal operation. A prerequisite for this is that the background lighting is switched off and that no peripheral devices are connected.

The device automatically switches to battery operation as soon as the mains supply is interrupted. When the mains supply is restored, the device automatically switches back to mains operation.

The battery symbol indicates the present charging level of the battery. 1 segment corresponds to approx. 25 % capacity. When the symbol flashes the battery must be charged (min. 4 hours). The charging period is extended if work is continued during charging. The battery is protected against overcharging.

The charging time of the storage battery amounts to approx. 6 hours. If the device continues to be operated during the charging process, the charging time is extended. The storage battery has a service life of approx. 1,000 charging/discharging cycles.

Note

The storage battery is also suitable for permanent mains operation.

→ In order to obtain the full nominal capacity we recommend that you discharge the storage battery at regular intervals (approx. every 4 weeks) through normal operation.

1.4 Disposal



In conformance with the European Directive 2002/96 EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of with domestic waste. This also applies to countries outside the EU, per their specific requirements.

→ Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

If the device is equipped with a storage battery:

The nickel metal hydride (NiMH) storage battery does not contain any heavy metals. However, it may not be disposed of with the normal refuse.

→ Observe the local regulations on the disposal of materials that are hazardous to the environment.

BBA422 / BBA425 / BBK422 **Operation**

2 Operation

2.1 Switching on and off

Switching on

→ Press ①.

The scale conducts a display test. When the weight display appears, the scale is ready to weigh.

Switching off

→ Press ①.

Before the display goes out, -OFF- appears briefly.

2.2 Zeroing / Zero point correction

Zeroing corrects the influence of slight changes on the load plate.

Manual

- 1. Unload scale.
- 2. Press *→*0*←*.

The zero display appears.

Automatic

In the case of scales that cannot be certified, the automatic zero point correction can be deactivated in the menu or the amount can be changed.

As standard, the zero point of the scale is automatically corrected when the scale is unloaded.

2.3 Simple weighing

- 1. Place weighing sample on scale.
- 2. Wait until the stability monitor **O** goes out.
- 3. Read weighing result.

OperationBBA422 / BBA425 / BBK422

2.4 Weighing with tare

2.4.1 Taring

→ Place the empty container on the scale and press 万丈.

The zero display and the symbol **NET** appear.

The tare weight remains saved until it is cleared.

2.4.2 Clearing the tare

→ Unload scale and press (Fig. 2).

The symbol **NET** goes out, the zero display appears.

or

→ Press (c).

The symbol **NET** goes out, the gross weight appears in the display.

If A . $\mathtt{CL-tr}$ is activated in the menu, the tare weight is automatically cleared as soon as the scale is unloaded.

2.4.3 Automatic taring

Prerequisite

A-tArE is activated in the menu under SCALE \rightarrow tArE, the symbol **T** flashes in the display.

The packaging material must be heavier than 9 display steps of the scale.

→ Place the container or packaging material on the scale.

The packaging weight is automatically saved as the tare weight, the zero display and the symbol **NET** appear.

BBA422 / BBA425 / BBK422

Operation

2.4.4 Chain tare

Prerequisite

The tare function CHAIn.tr is activated in the menu.

With this function it is possible to tare several times if, for example, cardboard is placed between individual layers in a container.

- Place the first container or packaging material on the scale and press <a>T.
 The packaging weight is automatically saved as the tare weight, the zero display and the symbol NET appear.
- 2. Weigh the weighing sample and read/print out the result.
- 3. Place the second container or packaging material on the scale and press or again.

The total weight on the scale is saved as the new tare weight. The zero display appears.

- 4. Weigh the weighing sample in the second container and read/print the result.
- 5. Repeat the last two steps for other containers.

2.5 Dynamic weighing

With dynamic weighing, the scale calculates the mean value from 56 weighing operations within 4 seconds.

With manual start Prerequisite

AVERAGE -> MAnual is selected in the menu.

The weighing sample must be heavier than 5 scale divisions.

dynamic result is then displayed with the symbol *.

- 1. Place the weighing sample on the scale and wait until it has stabilized.
- Press to start dynamic weighing.
 During dynamic weighing, horizontal segments appear in the display, and the
- 3. Unload the scale to be able to start a new dynamic weighing operation.

Operation BBA422 / BBA425 / BBK422

With automatic start

Prerequisite

AVERAGE -> AUTO is selected in the menu.

The weighing sample must be heavier than 5 scale divisions.

1. Place the weighing sample on the scale.

The scale starts the dynamic weighing automatically.

During dynamic weighing, horizontal segments appear in the display, and the dynamic result is then displayed with the symbol *.

2. Unload the scale to be able to perform a new dynamic weighing operation.

2.6 **Printing results**

If a printer or computer is connected to the scale, the weighing results can be printed out or sent to a computer.

→ Press (□→).

The display contents are printed out and transferred to the computer.

2.7 Cleaning



CAUTION!

Electric shock hazard!

▲ Before cleaning with a damp cloth, pull out the mains plug to disconnect the unit from the power supply.



CAUTION!

When the weighing pan has been removed, never clean the area under the load plate holder with a solid object!

This could damage the weighing cell.

Other cleaning information:

- · Use damp cloths.
- Do not use any acids, alkalis or strong solvents.
- Do not clean using a high-pressure cleaning unit or under running water.
- If very dirty, remove the weighing pan, protective cover (if present) and adjustable feet and clean these items separately.
- Follow all the relevant instructions regarding cleaning intervals and permissible cleaning agents.

BBA422 / BBA425 / BBK422

Settings in the menu

3 Settings in the menu

Settings can be changed and functions can be activated in the menu. This enables adaptation to individual weighing requirements.

The menu consists of 6 main blocks containing various submenus on several levels.

3.1 Operating the menu

3.1.1 Calling up the menu and entering the password

The menu differentiates between 2 operating levels: Operator and Supervisor. The Supervisor level can be protected by a password. When the device is delivered, both levels are accessible without a password.

Operator menu

- 1. Press (=>) and keep it pressed until COdE appears.
- 2. Press (again.

The menu item terminu appears. Only the submenu device is accessible.

Supervisor menu

- 1. Press (and keep it pressed until COdE appears.
- 2. Enter the password and confirm with .

 The first menu item SCALE appears.

Note No supervisor password has been defined when the device is first delivered. Therefore respond to the password inquiry with when you call up the menu for the first time. If a password has still not been entered after a few seconds, the scale returns to weighing mode.


Emergency password for Supervisor access to the menu

If a password has been issued for Supervisor access to the menu and you have forgotten it, you can still enter the menu:

ightharpoonup Press ightharpoonup 3 times and confirm with ightharpoonup.

Settings in the menu BBA422 / BBA425 / BBK422

3.1.2 Selecting and setting parameters



Scrolling on one level

- → Scroll forward: Press 🥱 .
- → Scroll back: Press (→0).

Activating menu items/ accepting selection

→ Press 🕞 .

Exiting menu

1. Press **①**.

The last menu item END appears.

- 2. Press 🕞.
 - The inquiry SAVE appears.
- 3. Confirm inquiry with to save the settings and return to weighing mode.
 -or-
- → Press (>T+) to discard changes and return to weighing mode.

BBA422 / BBA425 / BBK422 Settings in the menu

3.2 Overview

Level 1	Level 2	Level 3	Level 4	Level 5	Level 6	Page
SCALE	CAL					
	dISPLAY UNIt1 g, kg , oz,		, lb, t	lb, t		
		UNIt2	g, kg, oz	, lb, t		
		rESOLU				
		UNt.rOLL	ON, OFF			
	tArE	A-tArE	ON, OFF			21
		ChAIn.tr	ON, OFF			
		A.CL-tr	ON, OFF , 9)d		7
	ZErO	AZM	OFF; 0.5	d; 1 d; 2 d	l; 5 d; 10 d	21
	rEStArt	ON/ OFF	-			21
	FILtEr	VibrAt	LOW, MEd , HIGH,		21	
		Process	UNIVEr, do	OSING		
		StAbILI	FASt, StA	ndrd, PrECI	SE	
	FACt	tEMP	OFF, 1K,	2K, 3K, 5K		22
	Min.WEiG	ON/OFF	ON, OFF			22
	rESEt	SUrE?			22	
APPLIC	AVErAGE	OFF, AUto, MAnuAL			22	
	rESEt	SUrE?				22
tERMINL	dEVICE	SLEEP	OFF, 1 min	n, 3 min, 5	min, 15 min,	23
		PWr OFF	OFF, 1 min	n, 3 min , 5	min, 15 min,	
		b.LIGHt	ON, OFF,	5 sec, 10 s	sec, 30 sec,	
	ACCESS	SUPErVI	I			23
	rESEt	SUrE?				23

Settings in the menu BBA422 / BBA425 / BBK422

Level 1	Level 2	Level 3	Level 4	Level 5	Level 6	Page
COMMUNI	COM 1/COM 2	MOdE	Print	l	1	24
			A.Print			7
			CONTINU			
			dIALOG			7
			CONt.OLd			
			dIAL.OLd			
			dt-b	GrOSS	ON, OFF	
				tArE	ON, OFF	
				nEt	ON, OFF	
			dt-G	Gross	ON, OFF	_ _ _
				tArE	ON, OFF	
				nEt	ON, OFF	
			COnt-Wt		·	_
			2nd.dISP			
			InSt.Prn			
		PriNtEr	Type ASCII, LAbEL		24	
			tEmPLat	StdArd, tEMPLt1, tEMPLt2		
			ASCi.Fmt	LINE.FMt	MULtI	
					SINGLE	
					FIXEd	
				LENGtH	1 100	
				SEPArAt	, ;	
				Add LF	0 9	
		PArAMEt	bAUd	300 38400		25
			PAritY		nonE, 7 odd, EVEN , 8 EVEN	
			H.SHAKE	NO, XONXO	FF , nEt 422,	
			NEt.Addr	0 31		1
			ChECSuM	ON, OFF		
			Vcc	ON, OFF		1
		rSt.COMx	SUrE?	<u> </u>		25

BBA422 / BBA425 / BBK422 Settings in the menu

Level 2	Level 3	Level 4	Level 5	Level 6	Page
OPTION	EtH.NEt	th.NEt IP.AddrS, SUbNEt, GATEWAY			25
	USb	USb tESt			25
	diGitAL	IN 0 3			25
		OUT 0 3	AbV.Min, Un	ndErLd,	
		SEt.Pt 1			
SEt.Pt 2					
dEF.PrN	tEmPLt1/ tEMPLt2	LINE 1 LINE 8	-	•	26
			StArLN, Crl	LF, F FEEd	
tESt SC	intErN/Ext	ErN			27
KboArd					
display					
SNr					
LiSt					
rESEt.AL	SUrE?				
	OPTION def.PrN test sc KboArd display snr List	OPTION ETH.NET USB diGitAL dEF.PrN tempLt1/ tEMPLt2 tESt SC intErN/Ext KboArd dISPLAY SNr LiSt	OPTION ETH.NET IP.Addrs, USb USb tEST diGitAL IN 0 3 OUT 0 3 SET.PT 1 SET.PT 2 dEF.PTN tEMPLt1/tINE 1 tEMPLt2 LINE 8 tEST SC intern/Extern KboArd dISPLAY SNr LiSt	OPTION EtH.NET IP.Addrs, SUbNET, GATE USB USB TEST diGitAL IN 0 3 OFF, ZErO, Print, CLEZ OUT 0 3 OFF, StAbLI AbV.Min, Un OVERLD, STZ SET.PT 1 SET.PT 2 def.PrN tempLt1/ LINE 1 Not.USED, H GROSS, tARI STARLN, CRI test SC intern/Extern KboArd display SNr List	OPTION EtH.NET IP.Addrs, SUBNET, GATEWAY USB USB TEST diGitAL IN 0 3 OFF, ZErO, tAre, Print, CLEAr, Unit OUT 0 3 OFF, StAbLE, bEL.Min, AbV.Min, UnderLd, OVErLd, StAr, SET.PT 1 SET.PT 2 def.PrN templt1/ LINE 1 Not.USEd, HEAder, GrOSS, tAre, net, StArLN, Crlf, F FEEd test SC intern/Extern KboArd display SNr List

Settings in the menu BBA422 / BBA425 / BBK422

3.3 Scale settings (SCALE)

3.3.1 CAL – calibration (adjustment)

This menu item is not available for certified scales without internal calibration weight.

Internal	For scales with an internal calibration weight:
	1. Unload scale.
	2. Activate menu item CAL with (E). The scale calibrates with the internal calibration weight. —Int CAL— appears in the display. After calibration is completed, —done— appears briefly in the display, and the scale automatically changes to the next point of the scale menu.
External	For scales without an internal calibration weight:
	1. Unload scale.
	2. Activate menu item CAL with (). The scale determines the zero point. -0 - appears in the display. The calibration weight to be placed on the scale then flashes in the display.
	3. If necessary, change the weight value displayed with 🖅.
	4. Place the calibration weight on the scale and confirm with .
	The scale calibrates with the calibration weight loaded. After calibration is completed, $-donE-$ appears briefly in the display, and the scale automatically changes to the next point of the scale menu.

3.3.2 DISPLAY – weighing unit and display accuracy

UNIt1	Select weighing unit 1: g, kg, oz, lb, t		
UNIt2	Select weighing unit 2: g, kg, oz, lb, t		
rESOLU	Select readability (resolution), model-dependent		
UNt.rOLL	When UNT.roll is switched on, the weight value can be displayed in all available units with .		
Notes	 In the case of certified scales individual sub-items of the display menu item may not be available or only to a limited extent, depending on the respective country. On dual-range/dual interval scales, resolutions marked with I<->I 1/2 are divided up into 2 weighing ranges/intervals, e.g. 2 x 3000 d. 		

BBA422 / BBA425 / BBK422 Settings in the menu

3.3.3 TARE – tare function

A-tArE	Switching on/off automatic taring
CHAIn.tr	Switching on/off chain tare
A.CL-tr	Switching on/off automatic clearing of the tare weight when the load is removed from scale
	Possible settings: OFF, ON, 9d

3.3.4 ZERO – automatic zero update

AZM	On certified scales, this menu item does not appear.
	Switching on/off automatic zero update and selecting zeroing range.
	Possible settings: OFF; 0.5 d; 1 d; 2 d; 5 d; 10 d

3.3.5 RESTART – automatic saving of zero point and tare value

ON/OFF	When the Restart function is activated, the last zero point and tare value are saved.
	After switching off / on or after a power interruption, the device continues to work with
	the saved zero point and tare value.

3.3.6 FILTER – adaptation to the ambient conditions and the weighing type

VIbrAt	Adaptation to the ambient conditions
LOW	 Very steady and stable environment. The scale works very quickly, but is very sensitive to external influences.
MEd	Normal environment. The scale operates at medium speed.
Restless environment. The scale works more slowly, but is insensitive to influences.	
Process	Adaptation to the weighing process
UNIVEr	Universal setting for all weighing samples and normal weighing goods
dosing	Dispensing liquid or powdery weighing samples
StAbILI	Adjusting the stability detection
FASt	The scale operates very fast.
StAndrd	The scale operates at medium speed.
PrECISE	The scale operates with the greatest possible reproducibility.
	The slower the scale works, the greater the reproducibility of the weighing results.

Settings in the menu BBA422 / BBA425 / BBK422

3.3.7 FACT – automatic temperature-dependent adjustment

This menu item appears only on scales with an internal calibration weight.

tEMP	Defining the temperature difference for automatic calibration
OFF	Switching off automatic calibration in the case of a temperature difference
1K/2K/3K/5K	Automatic calibration in the case of a temperature change of 1 K, 2 K, 3 K or 5 K since the last adjustment

3.3.8 MIN.WEIG - minimum weight

This menu item appears only if the service technician has saved a minimum weight.

ON/OFF	Switching minimum weight function on/off
	If the weight on the scale falls below the stored minimum weight, an * appears on
	the display in front of the weight indicator.

3.3.9 RESET – resetting scale settings to factory settings

SUrE?	Confirmation inquiry
	Reset the scale settings to factory settings with
	Do not reset scale settings with T

3.4 Application settings (APPLICATION)

3.4.1 AVERAGE – determining the average weight for an unstable load

OFF	Calculating average weight switched off
AUtO	Calculating average weight with automatic start of the weighing cycle
MAnuAL	Calculating average weight with manual start of the weighing cycle via

3.4.2 RESET – resetting application settings to factory settings

SUrE?	Confirmation inquiry
	Reset the application settings to factory settings with
	Do not reset the application settings with

BBA422 / BBA425 / BBK422 Settings in the menu

3.5 Terminal settings (TERMINAL)

3.5.1 DEVICE – Sleep mode, energy-saving mode and display backlighting

SLEEP	This menu item only appears on devices in mains operation.
	When SLEEP is activated, the scale switches off display and backlighting after the time period set when not in use. The display and backlighting are switched on again at the press of a key or if the weight changes.
	Possible settings: OFF, 1 min, 3 min, 5 min
PWr OFF	This menu item only appears on devices in battery operation.
OFF/1 min/	When PWr OFF is activated, the device switches itself off automatically after approx. 3 minutes when not in use. Afterwards it has to be switched on using \bigcirc .
	Possible settings: OFF (switched off), 1 min, 3 min, 5 min, 15 min, 30 min
b.LIGHt	Switching the display backlighting on/off.
OFF/5 sec/	Setting whether and after which time the background lighting is to be switched off.
	Scales with a storage battery switch the background lighting off automatically by default when no action takes place at the scale for approx. 5 seconds.
	Possible settings: OFF (switched off), 5 sec, 10 sec, 30 sec, 1 min, ON (switched on)
Note	This menu item is accessible without a Supervisor password.

3.5.2 ACCESS – password for Supervisor menu access

SUPErVI	Password entry for Supervisor menu access
ENTER.C	Request to enter password
	→ Enter the password and confirm with 🕞
rEtYPE.C	Request to repeat the password entry
	→ Enter the password again and confirm with (=>
Notes	The password can consist of up to 4 characters.
	The key must not be part of the password. It is required for confirming the password.
	The key
	If you enter an impermissible code or make a typing error in the repetition, COdE.Err. appears in the display.

3.5.3 RESET – resetting terminal settings to the factory settings

SUrE?	Confirmation inquiry
	 Reset terminal settings to the factory settings with Do not reset the terminal settings with <a>T

Settings in the menu BBA422 / BBA425 / BBK422

3.6 Configuring interfaces (COMMUNICATION)

3.6.1 COM1/COM2 -> MODE – operating mode of the serial interface

Print	Manual data output to the printer with 🕞
A.Print	Automatic output of stable results to the printer (e.g. for series weighing operations)
CONTINU	Ongoing output of all weight values via the interface
dIALOG	Bi-directional communication via MT-SICS commands, control of the scale via PC
CONt.OLd	As per CONTINU, see above, but with 2 fixed blanks in front of the unit (compatible with Spider 1/2/3)
dIAL.OLd	As per dIALOG, see above, but with 2 fixed blanks in front of the unit (compatible with Spider 1/2/3)
dt-b	DigiTOL-compatible format.
GROSS	Transfer of the gross weight, identified with "B"
tArE	Transfer of the tare weight
nEt	Transfer of the net weight
dt-G	As per dt-b, see above, gross weight identified with "G"
COnt-Wt	TOLEDO Continuous mode
2nd.dISP	For connecting a second display (automatically activates the 5-V voltage supply at Pin 9)
InSt.Prn	Immediate manual data output to the printer with (not certifiable)

3.6.2 COM1/COM2 -> PRINTER – settings for protocol printout

This menu item only appears if the mode "Print" or "A.Print" is selected.

tYPE	Select the printer type
ASCII	ASCII printer, e.g. Sprinter 1
LabEL	Label printer, capable of printing graphics
tEmPLat	Selecting protocol printout
StdArd	Standard printout
tEmPLt1	Printout in accordance with Template 1
tEmPLt2	Printout in accordance with Template 2
ASCi.Fmt	Selecting formats for the protocol printout
LINE.Fmt	Line format: MULtI (multi-line), SINGLE (single-line) or FIXEd
LENGtH	Line length: 0 100 characters, appears only with line format MULtI or FIXEd
SEPArAt	Separator: , ; . / \ _ and space; appears only with line format SINGLE
Add LF	• Line feed: 0 9

BBA422 / BBA425 / BBK422 Settings in the menu

3.6.3 COM1/COM2 -> PARAMET – communication parameter

bAUd	Selecting baud rate: 300, 600, 1200, 2400, 4800, 9600, 19200, 38400 baud
PAritY	Selecting parity: 7 none, 8 none, 7 odd, 8 odd, 7 even, 8 even
H.SHAKE	Select handshake: NO, XONXOFF, NET 422 (network operation via the optional RS422/RS485 interface via 4-wire bus, only for COM1), NET 485 (network operation via the optional RS422/RS485 interface via 2-wire bus, only for COM1)
NET.Addr	Assigning network address: 0 31, only for NET 485
ChECSuM	Activating checksum byte (appears only in TOLEDO Continuous mode)
Vcc	Switching 5V voltage, e.g. for a bar code reader, on / off

3.6.4 COM1/COM2 -> RESET COM1/RESET COM2 - resetting serial interface to factory settings

SUrE?	Confirmation inquiry
	 Reset interface settings to factory settings with Do not reset the interface settings with

3.6.5 OPTION – configuring options

If no option is installed or is not yet configured, N.A. appears in the display.

EtH.NEt	Configuration of the Ethernet interface
IP.AddrS	Enter IP address
SUBNEt	Enter Subnet address
GATEWAY	Enter Gateway address
USb	Configuration of the USB interface
USb TEST	Test of the USB interface. After the test has been passed, rEAdY appears in the display.
diGitAL	Configuration of the digital inputs/outputs
IN 0 3	Configuring inputs 0 3
OFF	Input not assigned
ZErO	• Key (>0<-)
tArE	Key ⋟T↔
PriNt	• Key 🕞
CLEAr	• Key C Clear
UNIt	Key Units

Settings in the menu BBA422 / BBA425 / BBK422

OUT 0 3	Configuring outputs 0 3
OFF	Output not assigned
StAbLE	Stable weight value
bEL.MIN	Minimum weight not reached
AbV.MIN	Minimum weight reached or exceeded
UNdErLd	Insufficient load
OVErLd	Overload
StAr	Changed/calculated value
bEL.SP1	Setpoint 1 not reached
AbV.SP1	Setpoint 1 reached or exceeded
bEL.SP2	Setpoint 2 not reached
AbV.SP2	Setpoint 2 reached or exceeded
SEt.Pt1	Enter value for setpoint 1
SEt.Pt2	Enter value for setpoint 2

${\bf 3.6.6} \qquad {\bf DEF.PRN-configuring\ templates}$

tEMPLt1/tEMPLt2	Selecting Template 1 or Template 2
LINE 1 8	Select line
NOt.USEd	Line not used
HEAdEr	• Line as header. The contents of the header must be defined via an interface command, see Section 4.1.
GROSS	Gross weight
tArE	Tare weight
nEt	Net weight
StARLN	Line with ***
CrLF	Line feed (blank line)
F FEEd	Page feed

BBA422 / BBA425 / BBK422 Settings in the menu

3.7 Diagnosis and printing out of the menu settings (DIAGNOS)

tESt SC	
Internal	Testing scale with internal calibration weight
	-Int CAL- appears in the display during the test.
	 After completion of the test, ideally *d=0.0g briefly appears in the display, after which the scale changes to the next menu item KboArd.
External	Testing scale with external calibration weight
	1. The scale checks the zero point0 - appears in the display. The test weight flashes in the display.
	2. If necessary, change the weight value displayed with 🖅.
	3. Put the calibration weight on the scale and confirm with .
	4. The scale checks the calibration weight put on them.
	5. After the test is completed, the deviation from the last calibration briefly appears in the display, ideally *d=0.0g, after which the scale changes to the next menu item KboArd.
KboArd	Keyboard test
PUSH 1 6	Press the keys
	Note
	You cannot abort the keyboard test!
	If you have selected the menu item KboArd, you must press all keys.
display	Display test: The scale displays all functioning segments
SNr	Display of the serial number
LiSt	Printout of a list of all menu settings
rESEt.AL	Resetting all menu settings to the factory settings
SUrE?	Confirmation inquiry
	Reset all menu settings to the factory settings with ()
	Do not reset the menu settings with

Interface description BBA422 / BBA425 / BBK422

4 Interface description

4.1 SICS interface commands

The compact scales BBA422 / BBA425 / BBK422 support the command set MT-SICS (METTLER TOLEDO **S**tandard **I**nterface **C**ommand **S**et). With SICS commands, it is possible to configure, query and operate the scales from a PC. SICS commands are divided up into various levels.

4.1.1 Available SICS commands

	Command	Meaning				
LEVEL 0	@	Reset the scale				
	Ю	Inquiry of all available SICS commands				
	11	Inquiry of SICS level and SICS versions				
	12	Inquiry of scale data				
	13	Inquiry of scale software version				
	14	Inquiry of serial number				
	16	Inquiry of weighing parameters				
	S	Send stable weight value				
	SI	Send weight value immediately				
	SIR	Send weight value repeatedly				
	Z	Zero the scale				
	ZI	Zero immediately				
LEVEL 1	D	Write text into display				
	DW	Weight display				
	K	Keyboard check				
	SR	Send and repeat stable weight value				
	T	Tare				
	TA	Tare value				
	TAC	Clear tare				
	TI	Tare immediately				

In the case of Levels O and 1, these are commands which, if implemented, will function identically with all METTLER TOLEDO scales or weighing terminals.

In addition there are also further interface commands which apply either to the entire product series or to the particular application level. This and further information on the MT-SICS command set may be found in the MT-SICS Manual (Order Number 22 011 459 or at www.mt.com) or be obtained by request from your METTLER TOLEDO customer service representative.

BBA422 / BBA425 / BBK422

4.1.2 Requirements for communication between scale and PC

- The scale must be connected to the RS232, RS485, USB or Ethernet interface of a PC with a suitable cable.
- The interface of the scale must be set to "Dialog" mode, see Section 3.6.1.
- A terminal program must be available on the PC, e.g. HyperTerminal.
- The communication parameters baud rate and parity must be set in the terminal program and on the scale to the same values, see Section 3.6.3.

4.1.3 Notes on network operation via the optional interface RS422/485

Up to 32 scales can be networked with the optional RS422/485 interface. In network operation, the scales must be addressed from the computer before commands can be sent and weighing results received.

Address	Hex	ASCII
0	0x30	0
1	0x31	1
2	0x32	2
9	0x39	9
10	Ox3A	:
11	0x3B	;
31	0x4F	0

Des	scription of the steps	Host	Direction	Scale
1.	Host addresses the scale, e.g. with the address 3A hex.	<esc> :</esc>	>	
2.	Host sends a SICS command, e.g. SI	SI <crlf></crlf>	>	
3.	The scale confirms receipt of the command and sends the address back		<	<esc>:</esc>
4.	The scale responds to the command and returns control of the bus to the host		<	S_S45.02_kg <crlf></crlf>

Interface description BBA422 / BBA425 / BBK422

4.2 TOLEDO Continuous mode

4.2.1 TOLEDO Continuous commands

In TOLEDO Continuous mode the scale supports the following input commands:

Command	Meaning
P	Printing out the current result
T	Taring of the scale
Z	Zero setting of the display
С	Deleting of the current value

4.2.2 Output format in TOLEDO Continuous mode

Weight values are always transferred in TOLEDO Continuous mode in the following format:

	Statu	S		Field	1					Field 2							
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
STX	SWA	SWB	SWC	MSD	_	_	_	_	LSD	MSD	_	_	_	_	LSD	CR	CHK
Field	1		Cont-\	Vt: 6 di	gits fo	or the	weight	value	that is	s transt	ferred	witho	ut con	nma a	ind uni	t	1
Field	2		Cont-\	Vt: 6 di	gits fo	or the	tare w	eight t	hat is t	transfei	rred w	/ithout	comr	na an	d unit		
STX			ASCII character 02 hex, character for "start of text"														
SWA	SWB,	SWC	Status	words	A, B,	C, see	e belov	V									
MSD			Most s	Most significant digit													
LSD			Least	Least significant digit													
CR			Carria	Carriage Return, ASCII character OD hex													
CHK				arriage Return, ASCII character OD hex hecksum (2-complement of the binary sum of the 7 lower bits of all the characters sent eforehand incl. STX and CR)							e 7 Io	ent					

BBA422 / BBA425 / BBK422 Interface description

Status wor	d A							
		Status b	oit					
Function	Selection	6	5	4	3	2	1	0
Decimal	X00	0	1			0	0	0
position	XO					0	0	1
	Х					0	1	0
	0.X					0	1	1
	0.0X					1	0	0
	0.00X					1	0	1
	0.000X					1	1	0
	0.0000X					1	1	1
Numerical	X1	1		0	1			
increment	X2			1	0			
	X5			1	1			

Status word B	
Function/Value	Bit
Gross/Net: Net = 1	0
Sign: Negative = 1	1
Overload/Underload = 1	2
Movement = 1	3
lb/kg: kg = 1	4
1	5
Power up = 1	6

Status word C						
Function						
kg/lb	Bit					
0	1	0	1	0		
0	0	1	1	1		
0	0	0	0	2		
Print requ	uest = 1	"	1	3		
Extended	4					
1	5					
Tare mar	6					

Event and error messagesBBA422 / BBA425 / BBK422

5 Event and error messages

Error	Cause	Remedy
Display Dark	Back lighting set too dark	→ Set back lighting (b.LIGHt) brighter
	No mains voltage	→ Check mains
	Unit switched off	→ Switch on unit
	Mains cable not plugged in	→ Plug in mains plug
	Brief fault	→ Switch device off and back on again
Insufficient load	Load plate not on the scale	→ Place load plate on the scale
L J	Weighing range not reached	→ Set to zero
Overload	Weighing range exceeded	→ Unload scale
r 7		→ Reduce preload
	Result not yet stable	→ If necessary adjust vibration adapter or weigh dynamically
00	Function not permissible	→ Unload scale and set to zero
ר הם ד	Zeroing not possible with over- load or insufficient load	→ Unload scale
r _ u u _ ı		
Err 5	No calibration	 → Unplug the mains plug then plug it back in; switch unit off and then back on in battery mode → Calibrate scale → Call METTLER TOLEDO Service
Err 17	Printout not yet ended	→ End printout→ Repeat required action
Err 18	Switching the weighing unit impermissible during dynamic weighing	→ End dynamic weighing→ Switch weighing unit
Err 53	EAROM checksum error	 → Unplug the mains plug then plug it back in; switch unit off and then back on in battery mode → Call METTLER TOLEDO Service

Error	Cause	Remedy
Weight display unstable	Restless installation location	→ Adjust vibration adapter
	Draft	→ Avoid drafts
	Restless weighing sample	→ Dynamic weighing
	Contact between weighing pan and/or weighing sample and surroundings	→ Remedy contact
	Mains fault	→ Check mains
Incorrect weight display	Incorrect zeroing	→ Unload scale, set to zero and repeat weighing operation
	Incorrect tare value	→ Clear tare
	Contact between weighing pan and/or weighing sample and surroundings	→ Remedy contact
	Scale tilted	→ Level scale

BBA422 / BBA425 / BBK422

6 Technical data and accessories

6.1 Technical data

6.1.1 Type key

The compact scales BBA422 / BBA425 / BBK422 are available with various capacities and platforms that can be seen from the complete type designation.

Example

BBK422 – **3 XS** compact scale with capacity **3 kg** and **extra-small platform**BBA422 – **6 SM** compact scale with capacity **6 kg** and **small platform**

BBA425 - **35 LA** compact scale with capacity **35 kg** and large platform

6.1.2 General data

BBA422 / BBA425 / BBK42	2			
Applications	Weighing			
	Dynamic weighing			
Settings	Resolution selectable			
	Weighing unit selectable: g, kg, oz, lb, t			
	Taring function: manual, automatic, chain tare			
	Automatic zero point correction when the scale is switched on and during operation			
	Filter for adapting to the ambient conditions (vibration adapter)			
	• Filter for adapting to the weighing type, e.g. dispensing (weighing process adapter)			
	 Switch-off function, sleep mode for mains-operated devices, energy-saving mode for battery operation 			
	Display lighting			
Accuracy class OIML/NTEP	BBA4 III			
	• BBK4 II			
Display	LCD (liquid crystal display), digits 16 mm high, with back lighting			
Keypad	Pressure point membrane keypad			
	Scratch-proof labeling			
Housing	Diecast aluminum housing; chromium nickel steel weighing pan			
	Dimensions, see Page 37			
Protection Class (IEC 529,	BBA422 / BBK422 IP43 (not with Ethernet interface)			
DIN 40050, EN60529)	• BBA425 IP65			

BBA422 / BBA425 / BBK422 Technical data and accessories

BBA422 / BBA425 / BB	K422				
Mains connection	the nominal voltage): • 230 V, 50 Hz, 70 mA • 240 V, 50 Hz, 70 mA • 120 V, 60 Hz, 90 mA • 100 V, 50/60 Hz, 90 mA For battery operation: • Connection via mains adapter: 90 -	 230 V, 50 Hz, 70 mA 240 V, 50 Hz, 70 mA 120 V, 60 Hz, 90 mA 100 V, 50/60 Hz, 90 mA 			
	Infeed on the unit: 24 V, 1.3 A				
Battery operation	If the voltage supply is interrupted, the operation	unit automatically switches over to battery			
Ambient conditions	• Use	Indoor use only			
	Altitude	up to 2000 m			
	Temperature range BBA4	–10 +40 °C / 14 104 °F			
	Temperature range BBK4	+10 +30 °C / 50 86 °F			
	Installation/overvoltage category	II			
	Pollution degree	2			
	Relative humidity	Maximum relative humidity 80 % for temperatures up to 31 °C / 88 °F, decreasing linearly to 50 % relative humidity at 40 °C / 104 °F			
Interfaces	1 RS232 interface integrated				
	1 other optional interface possible				

Technical data and accessories BBA422 / BBA425 / BBK422

6.1.3 Weighing ranges and readability BBA4..

The compact scales BBA4.. with strain gauge weighing cells are supplied in the configuration $2 \times 3000 \, d$. Higher legibilities are available from the factory with the optional "Premium" weighing cells.


Capacity	Configuration			
	2 x 3000 d (sto	ındard)	1 x 6000 d (with optional "Premium" weighing cells)	
	Weighing ranges	Readability (certified)	Weighing range	Readability (certified)
3 kg	1.5 kg / 3 kg	0.5 g / 1 g	3 kg	0.5 g
6 kg	3 kg / 6 kg	1 g / 2 g	6 kg	1 g
15 kg	6 kg / 15 kg	2 g / 5 g	15 kg	2 g
35 kg	15 kg / 35 kg	5 g / 10 g	35 kg	5 g
60 kg	30 kg / 60 kg	10 g / 20 g	60 kg	10 g

6.1.4 Weighing ranges and readability BBK4..

Configurations up to 1 x 60.000 e are possible with the compact scales BBK4.. with MonoBloc technology. Certified compact scales BBK4.. are supplied as standard with an internal calibration weight.

Model	Weighing interval(s)	Readability d	Verification value e
BBK4 – 3 DXS	600 g / 3,100 g	0.01 g / 0.1 g	0.1 g
BBK4 – 3 XS	3,100 g	0.01	0.1 g
BBK4 – 6 DXS	1,200 g / 6,100 g	0.01 g / 0.1 g	0.1 g
BBK4 – 6 XS	6,100 g	0.01	0.1 g
BBK4 6 DSM	1,200 g / 6,100 g	0.1 g / 1 g	1 g
BBK4 – 6 SM	6,100 g	0.2 g	0.2 g
BBK4 – 15 DLA	3,500 g / 15,100 g	0.1 g / 1 g	1 g
BBK4 – 15 LA	15,100 g	0.5 g	0.5 g
BBK4 – 35 DLA	7,000 g / 35,100 g	0.1 g / 1 g	1 g
BBK4 – 35 LA	35,100 g	0.1 g	1 g

6.1.5 Dimensions



	A	В	С	D	E	F	G	Н	ı	K	L
XS ¹⁾	335	265	100	240	200	46	276	208	216	165	165
SM ¹⁾	335	265	100	240	200	46	276	208	216	_	_
LA ¹⁾	370	360	115	350	240	52	310	304	310	_	_

¹⁾ dimensions in mm

6.1.6 Net weights

Model	without battery	with battery	with internal calibration weight (without battery)
BBA422 SM	4.6 kg	5.3 kg	_
BBA425 – SM	4.7 kg	5.4 kg	_
BBA422 – LA	8.2 kg	8.9 kg	_
BBA425 – LA	8.3 kg	9.0 kg	_
BBK4 – XS	4.9 kg	5.6 kg	5.4 kg
BBK4 – SM	4.7 kg	5.4 kg	5.2 kg
BBK4 – LA	10.5 kg	11.2 kg	11.7 kg

Technical data and accessories BBA422 / BBA425 / BBK422

6.1.7 Interface connections

The compact scale can be fitted with a maximum of 2 interfaces. The following combinations are possible:

COM1	COM2	Note
RS232	_	
RS232	RS232	
RS485	RS232	COM1 can be optionally operated as RS422 or RS485
RS232	Ethernet	10BaseT, RJ45 (not for BBA425)
RS232	USB	USB 1.1, Type B
RS232	Digital I/O	4 x in, 4 x out, D-Sub 9

6.1.8 Assignment of the interface connections

Pin	RS232	RS422	RS485	Digital I/O
	(COM1/COM2)	(4-wire, COM1)	(2-wire, COM1)	(COM2)
1	_	_	_	GND
2	TxD1/2	TxD1-	TxD1-/RxD1-	OUT0
3	RxD1/2	RxD1-	_	OUT1
4	_	_	_	OUT2
5	GND	GND	GND	OUT3
6	_	_	_	INO
7	_	TxD1+	TxD1+/RxD1+	IN1
8	_	RxD1+	_	IN2
9	VCC	VCC	VCC	IN3

6.2 Accessories

Designation	Order number
Protective cover for small model	21 203 207
Protective cover for large model	21 203 206
Second display RS-PD/PASM	21 302 875
Second display ADI412	22 013 978
Second display ADI412-B, with backlighting	22 013 977
Relay box 4 for connection to digital I/O interface	22 011 967
Connection cable for relay box 4, length approx. 1.5 m	21 254 225
Printer Sprinter 1 Euro version	21 253 399
Printer Sprinter 1 UK version	21 253 745
Anti-theft device	00 229 175
RS232 cable for printer Sprinter 1, 1.8 m long	21 253 677
RS232 cable for PC, 1.8 m long	00 410 024
Adapter plate for glass windshield for BBK4XS	22 011 304
Special load plate 186 x 186 mm for BBK4XS	22 013 451

Appendix BBA422 / BBA425 / BBK422

7 Appendix

7.1 Information for certified scales in EC countries



Weighing instruments verified at the place of manufacture bear the preceding mark on the packing label and green "M" sticker on the descriptive plate. They may be set to work immediately.



Weighing instruments which are verified in two steps have no green "M" on the descriptive plate and bear the preceding identifaction mark on the packing label. The second step of the verification must be carried out by the approved METTLER TOLEDO Service or by the W&M authorities. Please contact METTLER TOLEDO Customer Service.

The first calibration step of the verification has been carried out at the manufacturing plant. It comprises all tests according to EN45501-8.2.2. Scales with analog connection to the weighing platform require an additional test according to EN45501-3.5.3.3. However, this test is not mandatory if the terminal bears the same serial number as the weighing platform.

If national regulations in individual countries limit the period of validity of the certification, the operator of such a scale is himself responsible for its timely re-certification.

7.2 Safety checks

The compact scales of the series BBA422 / BBA425 / BBK422 have been checked by accredited testing institutions. They have passed the safety checks listed below and carry the relevant test symbols. Production is subject to production monitoring by the inspection offices.

Country	Test symbol	Standard
Canada		CAN/CSA-C22.2 No. 1010.1-92
USA	c SP us	UL Std. No. 61010A-1
Other countries	CB Scheme	IEC/EN61010-1:2001
	(no identification)	

BBA422 / BBA425 / BBK422 Appendix

7.3 Table of Geo Values

For weighing instruments verified at the manufacturer's, the geo value indicates the country or geographical zone for which the instrument is verified. The geo value set in the instrument (e.g. "Geo 18") appears briefly after switch-on or is specified on a label.

Table **GEO VALUES 3000e** shows the geo values for European countries.

Table **GEO VALUES 6000e/7500e** shows the geo values for different gravitation zones.

7.3.1 GEO VALUES 3000e, OIML Class III (European Countries)

Geographical latitude	Geo value	Country	
46°22' – 49°01'	18	Austria	
49°30' – 51°30'	21	Belgium	
41°41' – 44°13'	16	Bulgaria	
42°24' – 46°32'	18	Croatia	
48°34' – 51°03'	20	Czechia	
54°34' – 57°45'	23	Denmark	
57°30' – 59°40'	24	Estonia	
59°48' – 64°00'	25*	Finland	
64°00' – 70°05'	26		
41°20' – 45°00'	17	France	
45°00' – 51°00'	19*		
47°00' – 55°00'	20	Germany	
34°48' – 41°45'	15	Greece	
45°45' – 48°35'	19	Hungary	
63°17' – 67°09'	26	Iceland	
51°05' – 55°05'	22	Ireland	
35°47' – 47°05'	17	Italy	
55°30' – 58°04'	23	Latvia	
47°03' – 47°14'	18	Liechtenstein	
53°54' – 56°24'	22	Lithuiania	
49°27' – 50°11'	20	Luxemburg	
50°46' – 53°32'	21	Netherlands	
57°57' – 64°00'	24*	Norway	
64°00' – 71°11'	26		
49°00' – 54°30'	21	Poland	
36°58' – 42°10'	15	Portugal	
43°37' – 48°15'	18	Romania	

Appendix BBA422 / BBA425 / BBK422

Geographical latitude	Geo value	Country	
47°44' – 49°46'	19	Slovakia	
45°26' – 46°35'	18	Slovenia	
36°00' – 43°47'	15	Spain	
55°20' – 62°00'	24*	Sweden	
62°00' – 69°04'	26		
45°49' – 47°49'	18	Switzerland	
35°51' – 42°06'	16	Turkey	
49°00' – 55°00'	21*	United Kingdom	
55°00' – 62°00'	23		

^{*} factory setting

BBA422 / BBA425 / BBK422 **Appendix**

7.3.2 GEO VALUES 6000e/7500e OIML Class III (Height \leq 1000 m)

Geograhical latitude	Geo value		
00°00' – 12°44'	5		
05°46' – 17°10'	6		
12°44' – 20°45'	7		
17°10' – 23°54'	8		
20°45' – 26°45'	9		
23°54' – 29°25'	10		
26°45' – 31°56'	11		
29°25' – 34°21'	12		
31°56' – 36°41'	13		
34°21' – 38°58'	14		
36°41' – 41°12'	15		
38°58' – 43°26'	16		
41°12' – 45°38'	17		
43°26' – 47°51'	18		
45°38' – 50°06'	19		
47°51' – 52°22'	20		
50°06' – 54°41'	21		
52°22' – 57°04'	22		
54°41' – 59°32'	23		
57°04' – 62°09'	24		
59°32' – 64°55'	25		
62°09' – 67°57'	26		
64°55' – 71°21'	27		
67°57' – 75°24'	28		
71°21' – 80°56'	29		
75°24' – 90°00'	30		

 Index
 BBA422 / BBA425 / BBK422

8 Index

A	0
Accessories	Operator menu15
Adjustment 20	Options5, 25
Alignment9	P
Ambient conditions	Password15
Applications 34	Power supply9
Average 13, 22	Protocol14
C	
Calibrate 20	R
Chain tare 13	Readability36
Continuous mode 30	Reset
D.	Application
Dimensions 27	Interface
Dimensions	Scale
Display goograpy 7	RS422/RS48529
Display accuracy	13422/1340329
Dynamic weighing 13, 22	S
E	Safety checks40
Error messages	Scales types5
F	Settings34
Filter 21	SICS commands28
1 IIIGI Z I	Supervisor menu15
I	Switching off
Interface protocol 30	Switching on11
Interfaces	Т
Configure24	Tare
Connections 38	Automatic12
K	Chain tare
Keyboard 8	Clear12
Reybodid	Terminal settings23
M	TOLEDO Continuous30
Mains connection35	Type key34
Menu	
Application22	W
Communication 24	Weighing ranges
Diagnosis27	Weighing unit20
Operation 15	Weight38
Overview 17	Z
Scale 20	Zeroing11
Terminal	Š
Menu structure	



22011381B

Subject to technical changes © Mettler-Toledo (Albstadt) GmbH 05/08 Printed in Germany 22011381B

Mettler-Toledo (Albstadt) GmbH

D-72458 Albstadt

Tel. ++49-7431-14 0, Fax ++49-7431-14 232

Internet: http://www.mt.com